

# Mobility client set up

## iPhone, iPad or Android

Android:	Go to Google Play Store
iPhone / iPad:	Go to App Store
Download the free software: (USA)	Fujifilm Synapse Mobility

## Set up the device

Server Name:	iuhsynmobility
Address:	iuhsynmobility.iuhealth.org:8443
Secure Connection:	ON

## Log on

Use your MHG username and password

---

## Internet browser access

Works with Mac/PC laptop or desktop

<https://iuhsynmobility.iuhealth.org:8443/pureweb/synapse/studybrowser>

Log on

Use your MHG username and password

---